



Home About Us Contact Us View Cart My Account FAQ

username

LOGIN

New Account »
Forgot Password?

Aluminum oxide



Advanced Search »

Ads by Google

R2S Microbial Air Sampler

Volumetric Cleanroom Air Samplers for detection of Bacteria and Mold
www.emtr2s.com

Particle Size Measurement

Concentrated system. No dilution. Zeta potential optional. Titration
www.Dispersion.com/Par

ISO 16232 Particle Sizing

Fully Automated Filter Scanning Saves Time - Accurate - Easy to Use
www.gt-vision.com

Identify Unknown Material

Powders, fibers, plastics, liquids, polymers, contaminations
www.EMSL.com

Chemistry Physical Chemistry

Sampling Efficiency Measurement Methods for Aerosol Samplers

Authors: [Jana S. Kesavan](#); [Jerold R. Bottiger](#); [Robert W. Doherty](#); [EDGEWOOD CHEMICAL BIOLOGICAL CENTER ABERDEEN PROVING GROUND MD RESEARCH AND TECHNOLOGY DIR](#)

Abstract: Aerosol samplers and concentrators are characterized at the U.S. Army Edgewood Chemical Biological Center. The characterization methods are summarized and include (1) monodisperse fluorescent/nonfluorescent Polystyrene Latex (PSL) microspheres with fluorometric analysis or Coulter Multisizer Analysis; (2) polydisperse solid **aluminum oxide** particles with Aerodynamic Particle Sizer (APS) analysis or Coulter Multisizer Analysis; (3) liquid fluorescent oleic acid particles with fluorometric analysis; and (4) bioparticles with Coulter Multisizer, culturing, enzyme linked immunosorbent assays (ELISA), polymerase chain reaction (PCR), and APS analyses. Samples are sent out for ELISA and PCR analysis. Aerosols are generated using several methods. The PSL microsphere aerosols are generated using a Collison nebulizer, sonic nozzle, Ink Jet Aerosol Generator (IJAG), and puffers.

Aluminum oxide particle aerosols are generated using a sonic nozzle. Fluorescent oleic acid particle aerosols are generated using a Vibrating Orifice Aerosol Generator (VOAG). Bioparticle aerosols are generated using puffers (metered dose inhaler devices), nebulizers, IJAG, and sonic nozzle.

Limitations: APPROVED FOR PUBLIC RELEASE

Description: Final rept. Jan 2000-Jun 2004

Pages: 21

Report Date: JAN 2005

Report Number: A529134



Keywords relating to this report:

- ✦ [AEROSOL GENERATORS](#)
- ✦ [AEROSOLS](#)
- ✦ [ALUMINUM OXIDES](#)
- ✦ [CHAIN REACTIONS](#)
- ✦ [EFFICIENCY](#)
- ✦ [FLUORESCENCE](#)
- ✦ [FLUOROMETERS](#)
- ✦ [IMMUNOASSAY](#)
- ✦ [MEASUREMENT](#)
- ✦ [MICROSPHERES](#)
- ✦ [NOZZLES](#)
- ✦ [OLEIC ACID](#)
- ✦ [ORIFICES](#)
- ✦ [PARTICLE SIZE](#)
- ✦ [POLYSTYRENE](#)
- ✦ [SAMPLERS](#)

Adobe PDF - \$13.95

Printed Format - \$23.95

ADD TO CART

Please check the box for the format you wish to order.

[Shipping Terms](#)
[About Electronic Delivery](#)

[Email This Abstract](#)

[« Back to search](#)

[Home](#) | [About Us](#) | [Contact Us](#) | [View Cart](#) | [Customer Service](#) | [Shipping Terms](#) | [Advanced Search](#) | [Privacy Policy](#) | [Restrictions on PDF Usage](#)

© 2001-2008 Storming Media LLC. All rights reserved.